MOBILITY SOLUTIONS

FOR YOUR SMALL & MEDIUM-SIZED MEDICAL DEVICES



Mediroll provides a variety of mobile solution to carry light weight medical devices or IT systems - Tablets, laptops - around the hospital campus.

Flat packed roll stand + Adapters to the devices you sell = A roll stand ready to go in 15 minutes!

Whether you are a distributor or a vendor, mediroll roll stands' common structure enables you to stock flat packed roll stands with a variety of interfaces and adapters to serve your entire line of mobile devices.





COST EFFECTIVE PRICE

Save money and hassle by ordering your roll stands from the same place



FLAT PACKAGING

Compact packaging make it easy to ship and store



QUICK ASSEMBLY

15 min to build from box to final product

CARTS ARE COMPATIBLE WITH THE DEVICES FROM THE FOLLOWING VENDORS

COMEN AIDAN **GIMA MIDMARK PHILIPS VIATOM** LIFESCOPE **PROGETTI SUNRAY CARDIOLINE VELA GE HEALTHCARE MINDRAY HOSPITECH SPACELABS ZUG** NIHON KOHDEN **FUKUDA FAZZINI WELCHALLYN LEPO DRAEGER SCHILLER MASIMO CRITICARE**

CONFIGURE YOUR OWN CART IN 5 STEPS

STEP 1

CHOOSE A BASE





- » Versatil, can be used while seating or standing
- » Supports heavier medical devices
- » Best suited for I.T. equipment





5 CASTERS BASE:


» Symetrical base allows for better maneuverability



STEP 2 CHOOSE COLUMN

4 CASTERS BASE:







5 CASTERS BASE:





STEP 3 CHOOSE CASTERS



COW-CST-SW-T75

STANDARD

4 CASTER BASE 2 with locks + 2 swivel



COW-CST-SW-T100

4" CASTERS

STABLE FOR HEAVY LOADS

2 with locks + 2 swivel



CST-75-DBL-DTL-SW


3" DOUBLE CASTERS

STANDARD 5 CASTER BASE

All 5 casters have locks


STEP 4 CHOOSE ACCESSORIES

Use the symbols bellow to guide you on the right accessories to choose for your cart



STEP 5 MEDICAL DEVICE ADAPTER

Specify the brand of your device, and we will provide you with a specific connector



CUSTOMIZATION

We can help you develop a custom solution to better suit your needs







5WRS-ACC

ACCESSORY RING

Necessary when fitting accessories to the 5 caster column







WASTE BASKET

ADJUSTABLE OXYGEN CYLINDER HOLDER

RECENT STUDY CASES









IT CARTS

STEP 1 **4 CASTER BASE**

STEP 2 CHOOSE COLUMN





STEP 3 **CHOOSE CASTERS**







STEP 4 CHOOSE ACCESSORIES

All I.T. carts needs to have the Work Surface

















5-STEP CART CONFIGURATION PROCESS

